


POM KEPITAL TX-21		Grade		
- KEPITAL TX-21		가	() Grade	.
-		,	가	Type
		ISO 1183	g/cm ³	1.39
		ISO 1133	g/10min	16
(23 , 50 %RH)		ISO 62	%	0.2
(1.8 MPa)		ISO 75		90
		UL 94	-	HB
		ISO 11357		165
		ISO 11359	X 10 ⁻⁵ /	13
		ISO 527	MPa	2,500
		ISO 527	MPa	57
		ISO 527	%	10
		ISO 527	%	33
		ISO 178	MPa	79
		ISO 178	MPa	2350
Charpy	(Notched, 23)	ISO 179/1eA	kJ/m ²	7.5
Charpy	(Notched, -30)	ISO 179/1eA	kJ/m ²	5.0
		IEC 60093	Ω	1X10 ¹⁶
		IEC 60093	Ω • cm	1X10 ¹⁴
		IEC 60243-1	kV/mm	-
(, t = 2 mm)		ISO 294-4	%	2.0
		ISO 1043	-	POM

Revision No : 7 (2020-04-03)



(: 0.1 %)

80 ~ 90 3 h ~ 4 h .

: 60 ~ 80 (140 ~ 176)
: 170 ~ 210 (338 ~ 410)

Mold	Bn(Nozzle)	B3(Metering)	B2(Compression)	B1(Feeding)	Hopper
60 ~ 80	180 ~ 210	190 ~ 200	180 ~ 190	170 ~ 180	60 ~ 80

: 150 mm/s ~ 200 mm/s
: Maximum 20 bar

2 30, T Tower 10 (04637) 94 15 (13901)
Tel: 02-6952-8377 Tel: 031- 436 -1300

1. 가 가 ()
Implants
2. 가 가